

## FITPACK® Container MK 2: DS-003



Left: Shows an opened FITPACK® Container MK 2, with a used syringe being placed inside.

Right: The Yellow FITPACK® MK 2 on the right has the code DS003Y the cut out shows both clean and used equipment.

The FITPACK® Container MK 2 DS003Y is the only personal sharps container that is certified by both the Australian Standard AS 4031) and the British Standard (BS 7320).



### Product Information

**Product Summary:** A compact high quality injection moulded container complete with integral moulded lid. Equipped with patented syringe disposal one-way flap syringe storage and disposal locking system.

**Description:** FITPACK® Container MK 2

**Product Code:** DS-003 (BLACK) OR DS-003Y (Yellow)

**Material:** Injection moulded polypropylene

**Nominal Capacity:** 250ml (safe recommended usable space)

**Capacity to spill line:** 330ml (Maximum container capacity)

**Container Height:** 153mm

**Container Length:** 80.5mm (Top) 78mm (Base)

**Container Width:** 38mm (Top) 36mm (Base)

**Container Shape:** Rectangular tube with rounded ends

**Syringe Quantity:**

1ml	2ml	5ml	10ml
10	3	NA	NA

**Quantity per carton:** 100


**Product Standards:** AS4031:1996, AS4939:2001, BS7320:1990

**Manufacturing Standards:** ISO 9001:2000.

**Application:** Injecting Drug Users and discrete remote use



Australian Standard  
AS 4031



BS 7320