

FITSTICK® Container: PA001/001



Left: Cut out showing used 1ml syringe inside FITSTICK® Container.

Right: FITSTICK® Container (Empty)



Product Information

Product Summary: A unique high quality injection moulded container has been designed to safely dispose of a single 1ml insulin syringe. It features superior impact, crush and needle penetration resistance. This product can only be used once as the container renders syringe unusable.

Description: FITSTICK® Container

Product Code: PA001/001

Material: Injection moulded Nylon

Nominal Capacity: Single 1ml insulin syringe

Container Length: 107mm

Dimensions: 26.5 x 14mm (top) 12 x 14.5mm (body).

Container Shape: Pen Shape

Quantity per carton: 1500

Product Standards: Manufactured to conform to AS4939:2001 for crush & AS4031: 1999 for needle penetration resistance.

Manufacturing Standards: ISO 9001:2000

Application: Diabetics and injecting drug users, First Aid and discrete remote use.